


Three sections:


I. Navigator:




Sections
Blue are defaulted open

Bands

Use the **Filter window** to search for a **specific item**, ie: respiratory rate, by using the drop down arrow.




The desired results will now appear on the flowsheet in **light green**.




Setting Preferences on Bands and Sections:


Bands have been pre-selected by **location**. **Collapse** and **expand** sections by a **right click** on the header, or click the -/+.



II. View Window:



Search for all items that have **specific characteristics**:



- Red=Critical Results
- Orange=High Results
- Blue=Low Results
- Brown=Abnormal Results
- Light Green=Results from the Filter Window
- Light Blue=Section Dividers
- Black=Final/Signed Results
- Purple=Documentation has not been Signed


Setting the Timescale:

Right click on the time header, select the desired timeframe.

Actual - results display within the timeframe specified.


Interval - results display in columns created for the frequency.

Insert Date/Time - customize the timeframe.



Note: the new column will disappear if the end user refreshes before documenting and signing in the column. The column will then have to be re-created.

III. Filter Window:




Changing the Date/Time Range:

Use the **arrows** on either side of the current blue date/time bar,




OR: **Right click** the date/time bar, select the desired timeframe.



Activity View:

Some **Tasks** have now been incorporated into Iview, ie: **Activity of Daily Living**. They may be accessed from the Task List, or launched from Iview using the Current Tasks & Overdue Tasks icons.

Note: if you close the patient's chart and then re-open it, Activity View will disappear. Click the task on the tasklist or the task icon on Iview to re-launch the task you were working on. Activity View cannot be customized. Contents are related to specific Task.




Using the Legend:

Click **Options** on the menu bar. Select **Show Legend**. The icons and their description display.



Review Results:



IView Buttons & Icons	
	Show/Hide Navigator
	Show/Hide Empty Columns
	Sign
	Cancel (prior to signing)
	Associate Monitor
	Actions (copy forward)
	Current Tasks
	Overdue Tasks
	Customize View
	Conditional Logic (more documentation is needed depending on choice made, ie: Pain Assessment tool used)
	Dynamic Groups (when you have multiples of the same things, ie: IVs)

Stony Brook University Medical Center



IView Overview for Student Nurses



Quick Reference Guide

See your Champion for Additional help using PowerChart

OR

Call the HELP desk At 4-HELP (444-4357)