


# Hill-Rom NaviCare Pillow Speaker



**LEGEND:**

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>1. Speaker</li> <li>2. Normal Nurse Call Button</li> <li>3. Pain Medication Call Button</li> <li>4. Bathroom Call Button</li> <li>5. MyStation Power Button</li> <li>6. Lights 1--Increase room brightness</li> <li>7. Lights 2—Decrease Room Brightness</li> <li>8. Rewind Button</li> <li>9. Fast Forward Button</li> <li>10. Mystation Navigation Button—UP</li> <li>11. Mystation Menu Navigation Button-- LEFT</li> <li>12. Mystation Menu Navigation Button—RIGHT</li> <li>13. Mystation Menu Navigation Button—BACK</li> </ul> | <ul style="list-style-type: none"> <li>14. MyStation Menu Navigation Button--DOWN</li> <li>15. MyStation Menu Navigation Button—SELECT</li> <li>16. Play/Pause Button</li> <li>17. Volume UP Button</li> <li>18. Volume DOWN Button</li> <li>19. Channel Up Button</li> <li>20. MyStation Menu Home Screen Button</li> <li>21. Channel Up Button</li> <li>22. TV Button</li> <li>23. Keypad</li> </ul> |
|--|--|